

**Justine Rhodes –Visium HD 3’
Clown Fish**

Samples received on 12/12/25

FF –Curls

Rhodes – 1 (4-5 curls)

Rhodes – 2 (4-5 Curls)

Rhodes – 3 (Section scraped from Slides 1,2,6,12)

Slides 4 – H& E stain

Moving Forward Adjustments Needed

Try to place curls for RIN inside the Eppendorf tube, not on the cap for the actual GE blocks, and make serial sections going well past the target regions.

1. Sections are placed on the wrong side of the slide
2. Location/Placing the sections are good
3. RIN quality needed is more than 7.0 ALL SAMPLES PASSED
4. Morphology the today's sections are much better and the processing just follow the same way
5. Without folding or curling and breakage.
6. We can discuss a time from now to perform the assay after the break.


RIN Assessment

Justin Rhodes (RhD)


12/12/2025	First Run		OCT Curls (Clown Fish)	
	S.No	Sample ID	RNA conc ng/μl	RIN
	1	RhD 1	5.16	9.4
	2	RhD 2	4.2	9.2
	3	RhD 3S	25.4	8.3
			Curls from slides	

DNA services



A1: RhD1
RQN: 9.4




A2: RhD2
RQN: 9.2




A3: RhD35
RQN: 8.3







H&E Stain



Slide 4

The red box within the images represents a 6.5×6.5 mm area, as required for Visium Spatial

H&E Stain



MF-2.1

The red box within the images represents a 6.5×6.5 mm area, as required for Visium Spatial